

Product Datasheet

Anti-Hepatitis B virus (HBcAg) V1822

Product Name	Anti-Hepatitis B virus (HBcAg)
Catalog Number	V1822
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive in ELISA with both adw and ayw Hb core antigens. Non-reactive with normal human serum
Specificity	Hepatitis B Virus core antigen
Matched Pairs	ELISA pairs (capture-conjugate) : V1822-V1823, V1822-V1822, V1822-V1824
Shelf Life	Four years from date of manufacture