

Product Datasheet

Anti-Hepatitis B virus (HBsAg) V1835

Product Name	Anti-Hepatitis B virus (HBsAg)
Catalog Number	V1835
Tested Applications	ELISA, WB, IHC
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with normal human serum
Specificity	Surface antigen; subtypes ad & ay (clone developed using native antigen)
Matched Pairs	ELISA pairs (capture-probe) : V1835-V1837 (capture-conjugate) : V1835-V1814
Shelf Life	Four years from date of manufacture