

Product Datasheet

Anti-Hepatitis C virus (NS4b) V1867

Product Name	Anti-Hepatitis C virus (NS4b)
Catalog Number	V1867
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with genotypes 1a and 1b
Specificity	NS4b, 1710-1730
Shelf Life	Four years from date of manufacture