

Product Datasheet

Anti-Hepatitis C virus (core) V1868

Product Name	Anti-Hepatitis C virus (core)
Catalog Number	V1868
Tested Applications	ELISA, WB, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with genotypes 1a, 1b, and 2a
Specificity	core, 33-38 aa
Matched Pairs	ELISA pairs (capture-conjugate): V1863-V1868, V1864-V1868
Shelf Life	Four years from date of manufacture