

Product Datasheet

Anti-Hepatitis C virus (NS5a) V1877

Product Name	Anti-Hepatitis C virus (NS5a)
Catalog Number	V1877
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with genotypes 1a and 1b. Non-reactive with genotype 2a
Specificity	NS5a
Shelf Life	Four years from date of manufacture