

Product Datasheet

Anti-HIV-1 (p24) V1941

Product Name	Anti-HIV-1 (p24)
Catalog Number	V1941
Tested Applications	ELISA, WB, IHC, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other HIV polypeptides
Specificity	HIV1 p24
Matched Pairs	ELISA pairs (capture-probe) : V1941-V1949
Shelf Life	Four years from date of manufacture