

Product Datasheet

Anti-HIV-1 (p24) V1949

Product Name	Anti-HIV-1 (p24)
Catalog Number	V1949
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	HIV1 p24 (native antigen)
Matched Pairs	ELISA pairs (capture-conjugate): V1946-V1949, V1941-V1949
Shelf Life	Four years from date of manufacture