

Product Datasheet

Anti-Human CD138 Monoclonal Antibody, Azide Free Clone B-A38 854.070.000

Product Name Anti-Human CD138 Monoclonal Antibody, Azide Free Clone B-A38

Catalog Number 854.070.000

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, IHC, FC

Immunogen U266 cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises Syndecan-1, the heparan sulfate proteoglycan, a 65-70 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells