

Product Datasheet

Anti-Ebola virus (NP) V2604

Product Name	Anti-Ebola virus (NP)
Catalog Number	V2604
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the NP of Sudan & Zaire sequences & VLPs in ELISA & WB
Specificity	NP of Sudan & Zaire sequences
Shelf Life	Four years from date of manufacture