

Product Datasheet

Anti-Ebola virus (VP40) V2612

Product Name	Anti-Ebola virus (VP40)
Catalog Number	V2612
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the VP 40 of Sudan & Zaire sequences & VLPs in ELISA & WB
Specificity	VP40 of Sudan & Zaire sequences
Matched Pairs	ELISA pairs (capture-conjugate) : V2614-V2612, V2613-V2612, V2615-V2612
Shelf Life	Four years from date of manufacture