

Product Datasheet

Anti-Entamoeba histolytica V3556

Product Name	Anti-Entamoeba histolytica
Catalog Number	V3556
Tested Applications	ELISA, LF IA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with E. histolytica & E. dispar in ELISA & LF IA
Specificity	E. histolytica & E. dispar
Matched Pairs	Self pairs
Shelf Life	Four years from date of manufacture