

Product Datasheet

Anti-Plasmodium falciparum (HRP II) V3569

Product Name	Anti-Plasmodium falciparum (HRP II)
Catalog Number	V3569
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Recognizes histidine rich protein II of P. falciparum
Matched Pairs	ELISA pairs (capture-conjugate) : V3569-V3564
Shelf Life	Four years from date of manufacture