

Product Datasheet

Anti-Human CD138 Monoclonal Antibody, PerCP Conjugated Clone B-A38 854.078.010

Product Name	Anti-Human CD138 Monoclonal Antibody, PerCP Conjugated Clone B-A38
Catalog Number	854.078.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PercP
Tested Applications	FC
Immunogen	U266 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Syndecan-1, the heparan sulfate proteoglycan, a 65-70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells