

Product Datasheet

Anti-Chikungunya virus V3585

Product Name	Anti-Chikungunya virus
Catalog Number	V3585
Tested Applications	ELISA, IFA (membrane staining)
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other Alphaviruses including EEE, WEE, VEE and other Arboviruses
Specificity	Chikungunya virus
Shelf Life	Four years from date of manufacture