

Product Datasheet

Anti-Japanese Encephalitis virus V3596

Product Name	Anti-Japanese Encephalitis virus
Catalog Number	V3596
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with TBE virus, Yellow Fever virus, West Nile virus, Dengue viruses (1-4), Zika virus and Chikungunya virus
Specificity	NS1 protein
Matched Pairs	ELISA pairs (capture-conjugate) : V3592-V3596, V3593-V3596
Shelf Life	Four years from date of manufacture