

Product Datasheet

Anti-Enterovirus 70 V3651

Product Name Anti-Enterovirus 70

Catalog Number V3651

Tested Applications ELISA, IFA

Brand ViroStat

Form/Appearance Purified from mouse ascites by protein A chromatography

Concentration ≥ 1 mg/ml

Storage 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

Isotype IgG3

Clonality Monoclonal

Purity $\geq 90\%$

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Non-reactive with numerous other Enteroviruses

Specificity Enterovirus type 70

Shelf Life Four years from date of manufacture