

## Product Datasheet

### Anti-Enterovirus V3654

<b>Product Name</b>	Anti-Enterovirus
<b>Catalog Number</b>	V3654
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Reactive with Rhinoviruses (9), Coxsackie B viruses (6), ECHO viruses (6), Coxsackie A viruses (4), EV71, and Poliovirus. Non-reactive with Influenza A, Influenza B, RSV, Adenovirus, Para 1, Para 2, Para 3, Chlamydia pneumoniae, Mycoplasma pneumoniae and Hepatitis A virus
<b>Specificity</b>	Pan-Enterovirus
<b>Matched Pairs</b>	Self pairs in ELISA & LF IA
<b>Shelf Life</b>	Four years from date of manufacture