

Product Datasheet

Anti-Gram positive bacteria V3811

Product Name	Anti-Gram positive bacteria
Catalog Number	V3811
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with <i>Listeria monocytogenes</i> (all serotypes), <i>Streptococcus pneumoniae</i> , <i>Staphylococcus aureus</i> , <i>Staphylococcus epidermidis</i> , <i>Enterococcus faecalis</i> , <i>Bacillus cereus</i> , <i>Bacillus subtilis</i> , Group A <i>Streptococcus</i>
Specificity	Lipoteichoic acid (LTA)
Matched Pairs	ELISA pairs (capture-conjugate): V3802-V3804, V3804-V3804, V3802-V3802,
Shelf Life	Four years from date of manufacture