

Product Datasheet

Anti-Coronavirus (NP) V3877

Product Name	Anti-Coronavirus (NP)
Catalog Number	V3877
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with NP of native SARS, MERS, by IFA & COVID-19 by ELISA. Non-reactive with Coronavirus strains 229E, NL63, & OC43
Specificity	NP of native SARS, MERS, & COVID-19
Shelf Life	Four years from date of manufacture