

Product Datasheet

Anti-Listeria monocytogenes V4212

Product Name	Anti-Listeria monocytogenes
Catalog Number	V4212
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with a membrane antigen of L. monocytogenes
Specificity	Membrane antigen of L. monocytogenes
Matched Pairs	ELISA pairs (capture-conjugate): V4211-V4212
Shelf Life	Four years from date of manufacture