

Product Datasheet

Anti-Listeria monocytogenes (p60) V4213

Product Name	Anti-Listeria monocytogenes (p60)
Catalog Number	V4213
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	L. monocytogenes secreted protein p60
Matched Pairs	ELISA pairs (capture-conjugate): V4213-V4214, V4214-V4213, V4215-V4213, V4213-V4216
Shelf Life	Four years from date of manufacture