

Product Datasheet

Anti-Metapneumovirus V4841

Product Name Anti-Metapneumovirus

Catalog Number V4841

Tested Applications ELISA, IFA

Alternative Names HMPV

Brand ViroStat

Form/Appearance Purified from mouse ascites by protein A chromatography

Concentration ≥ 1 mg/ml

Storage 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

Isotype IgG1

Clonality Monoclonal

Purity $\geq 90\%$

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Non-reactive with Respiratory Syncytial Virus, Influenza A, Influenza B, Adenovirus, & Parainfluenza viruses 1, 2, & 3

Specificity Matrix protein of human metapneumovirus

Shelf Life Four years from date of manufacture