

Product Datasheet

Anti-Norovirus (GI) V4924

| | |
|----------------------------|---|
| Product Name | Anti-Norovirus (GI) |
| Catalog Number | V4924 |
| Tested Applications | ELISA, WB |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Non-reactive with GII |
| Specificity | GI capsid protein |
| Matched Pairs | ELISA pairs (capture-conjugate) : V4921-V4924, V4922-V4924, V4923-V4924 |
| Shelf Life | Four years from date of manufacture |