

Product Datasheet

Anti-Norovirus (GII) V4936

Product Name	Anti-Norovirus (GII)
Catalog Number	V4936
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with GII.1 and GII.17 VLPS
Specificity	GII.1 and GII.17 VLPS
Shelf Life	Four years from date of manufacture