

Product Datasheet

Anti-Hantavirus (N) V4957

Product Name Anti-Hantavirus (N)

Catalog Number V4957

Tested Applications ELISA, IFA

Brand ViroStat

Form/Appearance Purified from mouse ascites by protein A chromatography

Concentration ≥ 1 mg/ml

Storage 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

Isotype IgG1

Clonality Monoclonal

Purity $\geq 90\%$

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Reactive with Hantaan, Dobrava, Seoul, & Saaremaa

Specificity Nucleocapsid protein (N) of several Hantaviruses tested

Shelf Life Four years from date of manufacture