

Product Datasheet

Anti-Astrovirus V4986

Product Name Anti-Astrovirus

Catalog Number V4986

Tested Applications ELISA

Brand ViroStat

Form/Appearance Purified from mouse ascites by protein A chromatography

Concentration ≥ 1 mg/ml

Storage 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.

Isotype IgG2a

Clonality Monoclonal

Purity $\geq 90\%$

Buffer 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide

Cross Reactivity Not determined

Specificity Capsid protein

Matched Pairs ELISA pairs (capture-conjugate) : V4985-V4987, V4988-V4987

Shelf Life Four years from date of manufacture