

Product Datasheet

Anti-Astrovirus V4988

Product Name	Anti-Astrovirus
Catalog Number	V4988
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Capsid protein
Matched Pairs	ELISA pairs (capture-conjugate) : V4988-V4986, V4988-V4987. Pairs in ELISA to detect type 1, other types not tested
Shelf Life	Four years from date of manufacture