

Product Datasheet

Anti-Bordetella pertussis V5111

Product Name	Anti-Bordetella pertussis
Catalog Number	V5111
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Dengue 1-4, TBE virus, West Nile virus, Yellow Fever virus, Japanese Encephalitis virus and Chikungunya virus
Specificity	NS1 protein
Matched Pairs	ELISA pairs (capture-conjugate) : V5122-V5123, V5122-V5124
Shelf Life	Four years from date of manufacture