

## Product Datasheet

### Anti-Zika virus (NS1) V5122

<b>Product Name</b>	Anti-Zika virus (NS1)
<b>Catalog Number</b>	V5122
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with Dengue 1-4, TBE virus, West Nile virus, Yellow Fever virus, Japanese Encephalitis virus and Chikungunya virus
<b>Specificity</b>	NS1 protein
<b>Matched Pairs</b>	ELISA pairs (capture-conjugate) : V5122-V5123, V5123-V5123, V5123-V5125, V5124-V5123
<b>Shelf Life</b>	Four years from date of manufacture