

Product Datasheet

Anti-Clostridium difficile (toxin A) V5147

Product Name	Anti-Clostridium difficile (toxin A)
Catalog Number	V5147
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with toxin B
Specificity	Toxin A of Clostridium difficile
Matched Pairs	ELISA pairs (capture-conjugate) : V5157-V5147, V5155-V5147
Shelf Life	Four years from date of manufacture