

Product Datasheet

Anti-Clostridium difficile (toxin B) V5156

Product Name	Anti-Clostridium difficile (toxin B)
Catalog Number	V5156
Tested Applications	LF, ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with other toxins
Specificity	Toxin B of Clostridium difficile
Matched Pairs	LF pairs: V5152-V5156
Shelf Life	Four years from date of manufacture