

Product Datasheet

Anti-Dengue virus (NS1) V5160

Product Name	Anti-Dengue virus (NS1)
Catalog Number	V5160
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	NS1 protein of all 4 serotypes
Matched Pairs	ELISA pairs (capture-conjugate) : V5162-V5163. V5162-V5164, V5162-V5165. LF pairs (capture-conjugate) : V5166-V5162
Shelf Life	Four years from date of manufacture