

Product Datasheet

Anti-M. tuberculosis (Hsp65/GroEI) V5178

Product Name	Anti-M. tuberculosis (Hsp65/GroEI)
Catalog Number	V5178
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with <i>M. bovis</i> , <i>M. avium</i> , <i>M. phlei</i> , <i>M. parafortuitum</i> , <i>Rhodococcus</i> sp., <i>B. subtilis</i> , <i>S. pneumoniae</i> , & <i>E. coli</i>
Specificity	Hsp65 (GroEI) of <i>M. tuberculosis</i>
Shelf Life	Four years from date of manufacture