

Product Datasheet

Anti-Dengue virus (Env III) V5184

Product Name	Anti-Dengue virus (Env III)
Catalog Number	V5184
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Domain III region of the envelope protein of Dengue virus serotypes 1-3
Shelf Life	Four years from date of manufacture