

Product Datasheet

Anti-Clostridium difficile (GDH) V5192

Product Name	Anti-Clostridium difficile (GDH)
Catalog Number	V5192
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with mammalian cells
Specificity	Glutamate dehydrogenase of C. difficile
Matched Pairs	ELISA pairs (capture-conjugate) : V5192-V5196
Shelf Life	Four years from date of manufacture