

Product Datasheet

Anti-h Papilloma virus L1 (HPV 6, 11, 16) V5306

| | |
|----------------------------|---|
| Product Name | Anti-h Papilloma virus L1 (HPV 6, 11, 16) |
| Catalog Number | V5306 |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Reactive with the L1 of types 6, 11, & 16 |
| Specificity | L1 of types 6, 11, & 16 |
| Shelf Life | Four years from date of manufacture |