

## Product Datasheet

### Anti-Mycoplasma genitalium V5601

<b>Product Name</b>	Anti-Mycoplasma genitalium
<b>Catalog Number</b>	V5601
<b>Tested Applications</b>	ELISA, IFA
<b>Brand</b>	ViroStat
<b>Form/Appearance</b>	Purified from mouse ascites by protein A chromatography
<b>Concentration</b>	≥ 1 mg/ml
<b>Storage</b>	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 90%
<b>Buffer</b>	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
<b>Cross Reactivity</b>	Non-reactive with Mycoplasma hominis, Mycoplasma pneumoniae, Ureaplasma urealyticum, Herpes virus type 2, Chlamydia trachomatis, Nesseria gonorrhoeae, and human cells
<b>Specificity</b>	M. genitalium
<b>Matched Pairs</b>	ELISA pairs (capture-conjugate) : V5601-V5602, V5602-V5601,
<b>Shelf Life</b>	Four years from date of manufacture