

# Product Datasheet

## Anti-Measles virus V6012

Product Name	Anti-Measles virus
Catalog Number	V6012
Tested Applications	ELISA, immunofluorescence
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Mumps virus, RSV, Para 1, Para 2, & Para 3
Specificity	Measles virus
Shelf Life	Four years from date of manufacture