

Product Datasheet

Anti-Neisseria gonorrhoeae V6142

Product Name	Anti-Neisseria gonorrhoeae
Catalog Number	V6142
Tested Applications	Indirect ELISA, IFA (whole cells)
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with 9 strains of <i>N. gonorrhoeae</i> , <i>N. lactamica</i> , & <i>N. meningitidis</i> . Non-reactive with <i>N. mucosa</i> , <i>N. perflava</i> , <i>N. sicca</i> , <i>Gardnerella vaginalis</i> , group B <i>Streptococcus</i> , <i>Candida albicans</i> , <i>Hemophilus influenzae</i> type B, <i>Chlamydia trachomatis</i> , <i>Trichomonas vaginalis</i> , <i>Mycoplasma genitalium</i> , <i>Ureaplasma urealyticum</i> & <i>Herpes simplex virus</i> type 2
Specificity	All strains of <i>N. gonorrhoeae</i> tested
Shelf Life	Four years from date of manufacture