

Product Datasheet

Anti-Neisseria gonorrhoeae V6145

| | |
|----------------------------|--|
| Product Name | Anti-Neisseria gonorrhoeae |
| Catalog Number | V6145 |
| Tested Applications | Indirect ELISA, IFA (whole cells) |
| Brand | ViroStat |
| Form/Appearance | Purified from mouse ascites by protein A chromatography |
| Concentration | ≥ 1 mg/ml |
| Storage | 2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 90% |
| Buffer | 0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide |
| Cross Reactivity | Reactive with 9 strains of <i>N. gonorrhoeae</i> , <i>N. lactamica</i> , & <i>N. meningitidis</i> . Non-reactive with <i>N. mucosa</i> , <i>N. perflava</i> , <i>N. sicca</i> , <i>Gardnerella vaginalis</i> , group B <i>Streptococcus</i> , <i>Candida albicans</i> , <i>Hemophilus influenzae</i> type B, <i>Chlamydia trachomatis</i> , <i>Trichomonas vaginalis</i> , <i>Mycoplasma genitalium</i> , <i>Ureaplasma urealyticum</i> & <i>Herpes simplex virus</i> type 2 |
| Specificity | All strains of <i>N. gonorrhoeae</i> tested |
| Shelf Life | Four years from date of manufacture |