

Product Datasheet

Anti-Salmonella V6346

Product Name	Anti-Salmonella
Catalog Number	V6346
Tested Applications	LF, ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	Salmonella paratyphi A. Target antigen unknown
Matched Pairs	Self-pairs in LF immunoassay
Shelf Life	Four years from date of manufacture