

Product Datasheet

Anti-Salmonella LPS core V6394

Product Name	Anti-Salmonella LPS core
Catalog Number	V6394
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with all salmonella O-serogroups tested: A, B, C1, C2, D, E1, E3, E4, F, G1, G2. Non-reactive with E. coli (mix), Klebsiella, Citrobacter, Pseudomonas, Yersinia, Shigella, Proteus & Legionella
Specificity	Common LPS core of all Salmonellae tested
Matched Pairs	ELISA pairs with LPS extraction (capture-conjugate) : V6391-V6394, V6392-V6394
Shelf Life	Four years from date of manufacture