

Product Datasheet

Anti-Mycoplasma hominis V6508

Product Name	Anti-Mycoplasma hominis
Catalog Number	V6508
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with M. genitalium, M. pneumoniae, and U. urealyticum
Specificity	M. hominis
Shelf Life	Four years from date of manufacture