

Product Datasheet

Anti-Mycoplasma pneumoniae (P1) V6535

Product Name	Anti-Mycoplasma pneumoniae (P1)
Catalog Number	V6535
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2b
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive with Ureaplasma urealyticum, Mycoplasma hominis, & Mycoplasma genitalium
Specificity	P1 adhesin
Matched Pairs	ELISA pairs (capture-conjugate): V6532-V6535, V6533-V6535, V6534-V6537, V6533-V6536
Shelf Life	Four years from date of manufacture