

# Product Datasheet

## Anti-Mycoplasma pneumoniae V6581

Product Name	Anti-Mycoplasma pneumoniae
Catalog Number	V6581
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with 70 & 80 kDa WB bands of M. pneumoniae. Non-reactive with uninfected cells
Specificity	70 & 80 kDa WB bands of M. pneumoniae
Shelf Life	Four years from date of manufacture