

Product Datasheet

Anti-Helicobacter pylori (OMP) V6621

Product Name	Anti-Helicobacter pylori (OMP)
Catalog Number	V6621
Tested Applications	ELISA, WB, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Non-reactive in ELISA with C. jejuni, E. coli (mix), Salmonella (mix), Shigella (mix), P. aeruginosa, Yersinia, and Citrobacter
Specificity	19kDa (OMP) of H. pylori
Shelf Life	Four years from date of manufacture