

Product Datasheet

Anti-Helicobacter pylori (flagellin) V6643

Product Name	Anti-Helicobacter pylori (flagellin)
Catalog Number	V6643
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the 54 kDa (flagellin) of H. pylori. Non-reactive in IFA with Giardia, Cryptosporidium, Enterococcus, E. Coli sp., Salmonella sp., Shigella sp., Campylobacter, & Vibrio
Specificity	54 kDa (flagellin) of H. pylori
Matched Pairs	Pairs for antigen detection (capture-conjugate) : V6641-V6643
Shelf Life	Four years from date of manufacture