

Product Datasheet

Anti-Helicobacter pylori (flagellin) V6644

Product Name	Anti-Helicobacter pylori (flagellin)
Catalog Number	V6644
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Reactive with the 54 kDa (flagellin) of <i>H. pylori</i> . Non-reactive in IFA with <i>Giardia</i> , <i>Cryptosporidium</i> , <i>Enterococcus</i> , <i>E. Coli</i> sp., <i>Salmonella</i> sp., <i>Shigella</i> sp., <i>Campylobacter</i> , & <i>Vibrio</i>
Specificity	54kDa (flagellin) of <i>H. pylori</i>
Matched Pairs	Pairs for antigen detection (capture-conjugate) : V6641-V6644
Shelf Life	Four years from date of manufacture