

Product Datasheet

Anti-Trichomonas vaginalis V7101

Product Name	Anti-Trichomonas vaginalis
Catalog Number	V7101
Tested Applications	ELISA, IFA (indirect)
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	T. vaginalis
Matched Pairs	ELISA pairs (capture-conjugate): V7181-V7101, V7102-V7101
Shelf Life	Four years from date of manufacture