

Product Datasheet

Anti-Trichomonas vaginalis V7161

Product Name	Anti-Trichomonas vaginalis
Catalog Number	V7161
Tested Applications	ELISA, IFA
Brand	ViroStat
Form/Appearance	Purified from mouse ascites by protein A chromatography
Concentration	≥ 1 mg/ml
Storage	2-8°C or -20°C. Avoid freeze-thaws. Maximum one freeze-thaw during storage.
Isotype	IgG3
Clonality	Monoclonal
Purity	≥ 90%
Buffer	0.01M Phosphate Buffered Saline, pH 7.2, containing 0.1% sodium azide
Cross Reactivity	Not determined
Specificity	T. vaginalis
Matched Pairs	ELISA pairs (capture-conjugate): V7141-V7161, V7171-V7161
Shelf Life	Four years from date of manufacture